

23andMe

HOME

REPORTS

TOOLS

RESEARCH

+

0

PP

play+47b3...

Overview

Scientific Details

Print

Cystic Fibrosis

Cystic fibrosis is a rare genetic disorder characterized by impaired lung and digestive function. A person must have two variants in the CFTR gene in order to have this condition.

play+47b3d3fa10, you have one of the variants we tested.

You could pass this variant on to your children.

1 variant detected

in the CFTR gene

1 variant not determined

We could not determine your result for one variant. This can be caused by random test error, other factors that interfere with the test, or if you have two copies of a tested variant.

How To Use This Test

This test does not diagnose any health conditions.

Please talk to a healthcare professional if this condition runs in your family, you think you might have this condition, or you have any concerns about your results.

[Review the Carrier Status tutorial](#)

[See Scientific Details](#)

Intended Uses

- Tests for **multiple variants** in the CFTR gene.
- To identify carrier status for cystic fibrosis.

Limitations

- Does **not test** for all possible variants for the condition.
- Does **not report** if someone has two copies of a tested variant.

Important Ethnicities

- This test is most relevant for people of **European, Hispanic/Latino**, and **Ashkenazi Jewish** descent.

You are a carrier.

You could pass this variant on to your children. People with only one variant are not expected to have cystic fibrosis.



We detected one variant for cystic fibrosis.

We could not determine your result for one of the tested variants.

Your results may be relevant for you if you're thinking about starting a family.

If you and your partner are both carriers, each child may have a **25% chance** of having this condition. Your relatives may also wish to consider testing if they plan to have children.



About Cystic Fibrosis

- When symptoms develop**

Symptoms typically develop during infancy.
- Typical signs and symptoms**

- Chronic cough
 - Lung infections
 - Pancreatic insufficiency
 - Malnutrition
 - Infertility in males


- Ethnicities most affected**

This condition is common in all ethnicities, but most common in people of European descent.
- How it's treated**

There is currently no known cure. Treatment focuses on managing symptoms and preventing complications such as lung infections and malnutrition.

Read more at: [Genetics Home Reference](#) [GeneReviews](#) [Mayo Clinic](#)

Consider talking to a healthcare professional if you are thinking about having children.



If you're starting a family, a genetic counselor can help you and your partner understand if additional testing might be appropriate.

Connect with a GC




Share your results with your family.

Share your report



If you have other concerns about your results, consult with a healthcare professional.

Print report



Learn more about this condition and connect with support groups.

Learn more

Cystic Fibrosis

Scientific Details

Cystic fibrosis is a rare genetic disorder characterized by impaired lung and digestive function. A person must have two variants in the CFTR gene in order to have this condition.


Cystic fibrosis is caused by variants in the CFTR gene.

CFTR

The CFTR [gene](#) contains instructions for making a [protein](#) called cystic fibrosis transmembrane conductance regulator. This protein helps control the salt and water balance of certain organs by allowing chloride ions to pass in and out of [cells](#). Certain [variants](#) in the CFTR gene disrupt this function, causing the lungs, pancreas, and other organs to produce abnormally thick mucus. This mucus can clog the respiratory tract, leading to signs and symptoms of cystic fibrosis.

Read more at [Genetics Home Reference](#)

Chromosome 7



You have one variant detected by this test. Your result for one of the tested variants could not be determined.

Variants Detected

[View All Tested Markers](#)

Marker Tested	Genotype*	Additional Information	
DeltaF508 Gene: CFTR Marker: i3000001	(-) Variant copy from one of your parents	 CTT Typical copy from your other parent	Biological explanation Typical vs. variant DNA sequence(s) Percent of 23andMe customers with variant References [5] ClinVar

*This test cannot distinguish which copy you received from which parent. This test also cannot determine whether multiple [variants](#), if detected, were inherited from only one parent or from both parents. This may impact how these variants are passed down.

23andMe always reports [genotypes](#) based on the 'positive' strand of the human [genome](#) reference sequence (build 37). Other sources sometimes report genotypes using the opposite strand.

Test Details

Indications for Use

The 23andMe PGS Carrier Status Test for Cystic Fibrosis is indicated for the detection of 28 [variants](#) in the CFTR [gene](#). This test is intended to be used to determine carrier status for cystic fibrosis in adults, but cannot determine if a person has two copies of a tested variant. The test is most relevant for people of [Ashkenazi Jewish](#), European, and Hispanic/Latino descent.

Special Considerations

- Symptoms of cystic fibrosis may vary depending on the variants involved.
- [ACMG](#) recommends carrier testing for cystic fibrosis for people of all ethnicities considering having children. This test includes 22 of 23 variants recommended for testing by ACMG.

Test Performance Summary

Carrier Detection Rate & Relevant Ethnicities

The "carrier detection rate" is an estimate of the percentage of carriers for this condition that would be identified by this test. Carrier detection rate differs by ethnicity and is provided only where sufficient data is available.

Ashkenazi Jewish	94%	[3, 5]
European	89%	[3, 5]
Hispanic/Latino	73%	[3, 5]
African American	65%	[3, 5]
Asian	55%	[3, 5]

Analytical Performance

Accuracy was determined by comparing results from this test with results from sequencing. Greater than 99% of test results were correct. While unlikely, this test may provide false positive or false negative results. For more details on the analytical performance of this test, refer to the package insert.

Warnings and Limitations

- This test does not cover all variants that could cause this condition.*
- This test does not diagnose any health conditions.
- Positive results in individuals whose ethnicities are not commonly associated with this condition may be incorrect. Individuals in this situation should consider genetic counseling and follow-up testing.
- Share results with your healthcare professional for any medical purposes.
- If you are concerned about your results, consult with a healthcare professional.

See the [Package Insert](#) for more details on use and performance of this test.

* Variants not included in this test may be very rare, may not be available on our genotyping platform, or may not pass our testing standards.

References

- [510\(k\): Illumina MiSeqDx Cystic Fibrosis Clinical Sequencing Assay](#)
- [American College of Obstetricians and Gynecologists Committee on Genetics. \(2011\). "ACOG Committee Opinion No. 486: Update on carrier screening for cystic fibrosis." Obstet Gynecol. 117\(4\):1028-31.](#)
- [Bobadilla JL et al. \(2002\). "Cystic fibrosis: a worldwide analysis of CFTR mutations--correlation with incidence data and application to screening." Hum Mutat. 19\(6\):575-606.](#)
- [Ong T et al. \(2001\). "Cystic Fibrosis and Congenital Absence of the Vas Deferens." \[Updated 2017 Feb 2\].](#)
- [Watson MS et al. \(2004\). "Cystic fibrosis population carrier screening: 2004 revision of American College of Medical Genetics mutation panel." Genet Med. 6\(5\):387-91.](#)

Change Log

Your report may occasionally be updated based on new information. This Change Log describes updates and revisions to this report.

Date

10/16/2015

Change

The variant 2789+5G>A (i4000320) was added to the report. Customers who have this variant will see this variant detected in their result, and see updated content in their report.

The carrier detection rate was updated for customers who self-report having Ashkenazi Jewish ancestry. The chances of still being a carrier were also updated for customers with no variants detected who self-report having Ashkenazi Jewish ancestry.

The chances of still being a carrier were updated for customers with no variants detected who self-report having European or African American ancestry.

10/16/2015

Due to improvements in data analysis, some customers who previously received a "Not Determined" result for one or more of the following genetic markers may see a genotype at these markers: i4000291, i4000294, i4000296, i4000300, i4000301, i4000305, i4000307, i4000308, i4000311, i4000314, i4000315, i4000317, i4000318, i4000321, i4000323, i4000324, i4000325. This may also update the overall report result for these customers.

Oct 21, 2015

Cystic Fibrosis report created.