


D-Bifunctional Protein Deficiency

DBPD is a rare genetic disorder. It is characterized by abnormal muscle tone, developmental disability, seizures, and early death. A person must have two variants in the HSD17B4 gene in order to have this condition.

Overview Scientific Details

Jamie, you **do not have the variants** we tested.

You could still have a variant not covered by this test.



How To Use This Test

This test does not diagnose any health conditions.

Please talk to a healthcare professional if this condition runs in your family, you think you might have this condition, or you have any concerns about your results.

[Review the Carrier Status tutorial](#)
[See Scientific Details](#)

+ Intended Uses

- Tests for **multiple variants** in the HSD17B4 gene.
- To identify carrier status for DBPD.

- Limitations

- Does **not test** for all possible variants for the condition.
- Does **not report** if someone has two copies of a tested variant.

🌐 Important Ethnicities

- This test does **not** include the majority of HSD17B4 variants that cause DBPD in any ethnicity.

You are likely not a carrier.


This result may be less relevant for you because the variants that cause DBPD are rarely found in people of your ethnicity.



We ruled out the tested variants for DBPD. These variants are very rare in all ethnicities.

You still have a chance of being a carrier for DBPD.

We cannot estimate your chances because this condition is rare and not well studied in your ethnicity.



About D-Bifunctional Protein Deficiency

Also known as: Peroxisomal Bifunctional Enzyme Deficiency. DBPD is a form of leukodystrophy.

📅 When symptoms develop

Symptoms typically develop at birth or during infancy.

🏠 Typical signs and symptoms

- Abnormal muscle tone
- Seizures
- Developmental disability
- Hearing and vision loss
- Distinctive facial features
- Early death

👥 Ethnicities most affected


This condition is rare in all ethnicities.

🏥 How it's treated

There is currently no known cure. Treatment focuses on managing symptoms and preventing complications.

Read more at: [Genetics Home Reference](#) [GeneReviews](#) [NIH Office of Rare Diseases Research](#)

Consider talking to a healthcare professional if you are concerned about your results.



If you're starting a family, a genetic counselor can help you and your partner understand if additional testing might be appropriate.

Connect with a GC



Share your results with a healthcare professional.

Print report



Learn more about this condition and connect with support groups.

Learn more



Give the gift of DNA discovery.

Gift a kit

Refer friends, earn rewards.

Get reward

ANCESTRY

- Ancestry Overview
- All Ancestry Reports
- Ancestry Composition
- DNA Relatives
- Order Your DNA Book

HEALTH & TRAITS

- Health & Traits Overview
- All Health & Traits Reports
- My Health Action Plan
- Health Predisposition
- Pharmacogenetics
- Carrier Status
- Wellness
- Traits

RESEARCH

- Research Overview
- Surveys and Studies
- Edit Answers
- Publications

FAMILY & FRIENDS

- View all DNA Relatives
- Family Tree
- Your Connections
- GrandTree
- Advanced DNA Comparison

D-Bifunctional Protein Deficiency

DBPD is a rare genetic disorder. It is characterized by abnormal muscle tone, developmental disability, seizures, and early death. A person must have two variants in the HSD17B4 gene in order to have this condition.

[Overview](#)
[Scientific Details](#)


DBPD is caused by variants in the HSD17B4 gene.

HSD17B4

The HSD17B4 gene contains instructions for making an enzyme called D-bifunctional protein. This enzyme is found in peroxisomes, where it helps break down molecules called fatty acids to make energy. Certain variants in HSD17B4 disrupt this function, causing a harmful buildup of fat molecules inside of cells.

Read more at [Genetics Home Reference](#)



Chromosome 5



You have no variants detected by this test.

Variants Detected

View All Tested Markers

Marker Tested	Your Genotype*	Additional Information
G16S Gene: HSD17B4 Marker: i5007145	G Typical copy from one of your parents 	G Typical copy from your other parent <ul style="list-style-type: none"> Biological explanation Typical vs. variant DNA sequence(s) Percent of 23andMe customers with variant References [2, 3, 5, 7] ClinVar
N457Y Gene: HSD17B4 Marker: i5007146	A Typical copy from one of your parents 	A Typical copy from your other parent <ul style="list-style-type: none"> Biological explanation Typical vs. variant DNA sequence(s) Percent of 23andMe customers with variant References [2, 4, 6] ClinVar

*This test cannot distinguish which copy you received from which parent. This test also cannot determine whether multiple variants, if detected, were inherited from only one parent or from both parents. This may impact how these variants are passed down.

23andMe always reports genotypes based on the 'positive' strand of the human genome reference sequence (build 37). Other sources sometimes report genotypes using the opposite strand.

Test Interpretation

Post-test carrier risk for DBPD is the chance of still being a carrier for the condition if you do not have the variants tested. This chance depends on how common it is to be a carrier for DBPD and whether the variants we tested tend to be found in people of your ethnicity.

Because you do not have the variants we tested, your chances of still being a carrier are lower than for someone who has not been tested. However, we cannot provide an exact estimate because the information needed to calculate post-test carrier risk is not available for your ethnicity.

Test Details

Indications for Use

The 23andMe PGS Carrier Status Test for D-Bifunctional Protein Deficiency is indicated for the detection of two variants in the HSD17B4 gene. This test is intended to be used to determine carrier status for DBPD in adults, but cannot determine if a person has two copies of a tested variant.

Special Considerations

- This test does not include the majority of HSD17B4 variants that cause DBPD in any ethnicity.
- There are currently no professional guidelines in the U.S. for carrier testing for this condition.

Test Performance Summary

Carrier Detection Rate & Relevant Ethnicities

The "carrier detection rate" is an estimate of the percentage of carriers for this condition that would be identified by this test. Carrier detection rate differs by ethnicity and is provided only where sufficient data is available.

Worldwide	35%	[2]
-----------	-----	-------

Analytical Performance

Accuracy was determined by comparing results from this test with results from sequencing. Greater than 99% of test results were correct. While unlikely, this test may provide false positive or false negative results. For more details on the analytical performance of this test, refer to the package insert.

Warnings and Limitations

- This test does not cover all variants that could cause this condition.*
- This test does not diagnose any health conditions.
- Positive results in individuals whose ethnicities are not commonly associated with this condition may be incorrect. Individuals in this situation should consider genetic counseling and follow-up testing.
- Share results with your healthcare professional for any medical purposes.
- If you are concerned about your results, consult with a healthcare professional.

See the [Package Insert](#) for more details on use and performance of this test.

* Variants not included in this test may be very rare, may not be available on our genotyping platform, or may not pass our testing standards.

References

- Ferdinandusse S et al. (2006). "Clinical and biochemical spectrum of D-bifunctional protein deficiency." *Ann Neurol.* 59(1):92-104. ^
- Ferdinandusse S et al. (2006). "Mutational spectrum of D-bifunctional protein deficiency and structure-based genotype-phenotype analysis." *Am J Hum Genet.* 78(1):112-24. ^
- Gin YM et al. (2000). "Human peroxisomal multifunctional enzyme type 2. Site-directed mutagenesis studies show the importance of two protic residues for 2-enoyl-CoA hydratase 2 activity." *J Biol Chem.* 275(7):4965-72. ^
- Tsuchida S et al. (2012). "Hydratase activities of green fluorescent protein tagged human multifunctional enzyme type 2 hydratase domain and its variants." *J Oleo Sci.* 61(8):443-50. ^
- Van Grunsven EG et al. (1998). "Peroxisomal D-hydroxyacyl-CoA dehydrogenase deficiency: resolution of the enzyme defect and its molecular basis in bifunctional protein deficiency." *Proc Natl Acad Sci U S A.* 95(5):2128-33. ^
- Van Grunsven EG et al. (1999). "Enoyl-CoA hydratase deficiency: identification of a new type of D-bifunctional protein deficiency." *Hum Mol Genet.* 8(8):1509-16. ^
- Van Grunsven EG et al. (1999). "Peroxisomal bifunctional protein deficiency revisited: resolution of its true enzymatic and molecular basis." *Am J Hum Genet.* 64(1):99-107. ^
- Vanderver A et al. (2014). "Leukodystrophy Overview." ^

Change Log

Your report may occasionally be updated based on new information. This Change Log describes updates and revisions to this report.

Date	Change
Feb. 18, 2016	Due to improvements in data analysis, some customers who previously received a "Not Determined" result for i5007146 may see a genotype at this marker. This may also update the overall report result for these customers.
Oct. 21, 2015	D-Bifunctional Protein Deficiency report created.



Give the gift of DNA discovery.

Gift a kit

Refer friends, earn rewards.

Get reward

ANCESTRY

Ancestry Overview
 All Ancestry Reports
 Ancestry Composition
 DNA Relatives
 Order Your DNA Book

HEALTH & TRAITS

Health & Traits Overview
 All Health & Traits Reports
 My Health Action Plan
 Health Predisposition
 Pharmacogenetics
 Carrier Status
 Wellness
 Traits

RESEARCH

Research Overview
 Surveys and Studies
 Edit Answers
 Publications

FAMILY & FRIENDS

View all DNA Relatives
 Family Tree
 Family Connections
 GrandTree
 Advanced DNA Comparison