


Limb-Girdle Muscular Dystrophy Type 2D

LGMD2D is a rare genetic disorder. It is characterized by muscle weakness that worsens over time as well as heart and lung problems. A person must have two variants in the SGCA gene in order to have this condition.

Overview Scientific Details

Jamie, you **do not have the variant we tested.**

You could still have a variant not covered by this test.



How To Use This Test

This test does not diagnose any health conditions.

Please talk to a healthcare professional if this condition runs in your family, you think you might have this condition, or you have any concerns about your results.

[Review the Carrier Status tutorial](#)
[See Scientific Details](#)

+ Intended Uses

- To test for the [R77C variant](#) in the [SGCA gene](#).
- To identify [carrier](#) status for LGMD2D.

- Limitations

- Does **not test** for all possible variants for the condition.
- Does **not report** if someone has two copies of a tested variant.
- Does **not cover** other types of LGMD.

🌐 Important Ethnicities

- This test is expected to identify the majority of carriers of **Finnish** descent.

You are likely not a carrier.


This result may be less relevant for you because the variants that cause LGMD2D are rarely found in people of your ethnicity.



We ruled out the tested variant for LGMD2D. This variant is very rare in all ethnicities.

You still have a chance of being a carrier for LGMD2D.

We cannot estimate your chances because this condition is rare and not well studied in your ethnicity.



About Limb-Girdle Muscular Dystrophy Type 2D

Also known as: Alpha-Sarcoglycanopathy

📅 When symptoms develop

Symptoms typically develop between early childhood and adolescence.

🌡️ Typical signs and symptoms

- Wasting of arm and leg muscles closest to the torso
- Large calf muscles
- Curvature of the spine
- Heart and lung problems
- Shortened lifespan

👥 Ethnicities most affected


This condition is rare in all ethnicities.

🩺 How it's treated

There is currently no known cure. Therapy focuses on maintaining muscle function, preventing skeletal problems, and monitoring heart and lung function.

Read more at: [Genetics Home Reference](#) [GeneReviews](#)

Consider talking to a healthcare professional if you are concerned about your results.



If you're starting a family, a genetic counselor can help you and your partner understand if additional testing might be appropriate.

[Connect with a GC](#)



Share your results with a healthcare professional.

[Print report](#)



Learn more about this condition and connect with support groups.

[Learn more](#)



Give the gift of DNA discovery.

[Gift a kit](#)

Refer friends, earn rewards.

[Get reward](#)

ANCESTRY

- Ancestry Overview
- All Ancestry Reports
- Ancestry Composition
- DNA Relatives
- Order Your DNA Book

HEALTH & TRAITS

- Health & Traits Overview
- All Health & Traits Reports
- My Health Action Plan
- Health Predisposition
- Pharmacogenetics
- Carrier Status
- Wellness
- Traits

RESEARCH

- Research Overview
- Surveys and Studies
- Edit Answers
- Publications

FAMILY & FRIENDS

- View all DNA Relatives
- Family Tree
- Your Connections
- GrandTree
- Advanced DNA Comparison

Limb-Girdle Muscular Dystrophy Type 2D

LGMD2D is a rare genetic disorder. It is characterized by muscle weakness that worsens over time as well as heart and lung problems. A person must have two variants in the SGCA gene in order to have this condition.

Overview **Scientific Details**


LGMD2D is caused by variants in the SGCA gene.

SGCA


The SGCA gene contains instructions for making one part of a group of [proteins](#). These proteins, called the sarcoglycan protein complex, are found in muscle tissue where they help strengthen and protect muscle fibers. Certain [variants](#) in the SGCA gene prevent the protein complex from working properly.

Read more at [Genetics Home Reference](#)*

Chromosome 17



You have no variants detected by this test.

Variants Detected		View All Tested Markers	
Marker Tested	Your Genotype*	Additional Information	
R77C Gene: SGCA Marker: rs28933693	C Typical copy from one of your parents	 C Typical copy from your other parent	<ul style="list-style-type: none"> Biological explanation Typical vs. variant DNA sequence(s) Percent of 23andMe customers with variant References [1, 2, 3, 4] ClinVar*

*This test cannot distinguish which copy you received from which parent. This test also cannot determine whether multiple [variants](#), if detected, were inherited from only one parent or from both parents. This may impact how these variants are passed down.

23andMe always reports [genotypes](#) based on the 'positive' strand of the human [genome](#) reference sequence (build 37). Other sources sometimes report genotypes using the opposite strand.

Test Interpretation

This report provides an estimate of the chances of still being a carrier for people who do not have the [variant\(s\)](#) tested. This is known as the **post-test carrier risk**.

Post-test carrier risk is based on the average chance of being a carrier for a given ethnicity and the carrier detection rate of the test for a given ethnicity.

[View technical article on estimating post-test carrier risk.](#)

Post-Test Carrier Risk

This report provides an estimate of the post-test carrier risk for people of Finnish descent only.

- For people of partial Finnish descent, post-test carrier risk is less than that for those who are fully Finnish. The exact post-test risk depends on how much Finnish ancestry a person has.
- Post-test risk for other ethnicities cannot be provided because sufficient data is not available.

Post-test carrier risk for relevant ethnicities

Finnish	1 in 5,500	[4]
---------	------------	-------

Test Details

Indications for Use

The 23andMe PGS Carrier Status Test for Limb-Girdle Muscular Dystrophy Type 2D is indicated for the detection of the [R77C variant](#) in the [SGCA gene](#). This test is intended to be used to determine carrier status for LGMD2D in adults, but cannot determine if a person has two copies of a tested variant.

Special Considerations

- Symptoms can vary greatly in people with this condition, and can be mild in some cases.
- There are currently no professional guidelines in the U.S. for carrier testing for this condition.

Test Performance Summary

Carrier Detection Rate & Relevant Ethnicities

The "carrier detection rate" is an estimate of the percentage of carriers for this condition that would be identified by this test. Carrier detection rate differs by ethnicity and is provided only where sufficient data is available.

Finnish	95%	[4]
---------	-----	-------

Analytical Performance

Accuracy was determined by comparing results from this test with results from sequencing. Greater than 99% of test results were correct. While unlikely, this test may provide false positive or false negative results. For more details on the analytical performance of this test, refer to the package insert.

Warnings and Limitations

- This test does not cover all variants that could cause this condition.*
- This test does not diagnose any health conditions.
- Positive results in individuals whose ethnicities are not commonly associated with this condition may be incorrect. Individuals in this situation should consider genetic counseling and follow-up testing.
- Share results with your healthcare professional for any medical purposes.
- If you are concerned about your results, consult with a healthcare professional.

See the [Package Insert](#) for more details on use and performance of this test.

* Variants not included in this test may be very rare, may not be available on our genotyping platform, or may not pass our testing standards.

References

- [Bartoli M et al. \(2008\). "Mannosidase I inhibition rescues the human alpha-sarcoglycan R77C recurrent mutation." Hum Mol Genet. 17\(9\):1214-21. ^](#)
- [Carrié A et al. \(1997\). "Mutational diversity and hot spots in the alpha-sarcoglycan gene in autosomal recessive muscular dystrophy \(LGMD2D\)." J Med Genet. 34\(6\):470-5. ^](#)
- [Draviam RA et al. \(2006\). "Alpha-sarcoglycan is recycled from the plasma membrane in the absence of sarcoglycan complex assembly." Traffic. 7\(7\):793-810. ^](#)
- [Hackman P et al. \(2005\). "Enrichment of the R77C alpha-sarcoglycan gene mutation in Finnish LGMD2D patients." Muscle Nerve. 31\(2\):199-204. ^](#)
- [Pegoraro E et al. \(2000\). "Limb-Girdle Muscular Dystrophy Overview." \[Updated 2012 Aug 30\]. ^](#)

Change Log

Your report may occasionally be updated based on new information. This Change Log describes updates and revisions to this report.

Date	Change
Oct. 21, 2015	Limb-Girdle Muscular Dystrophy Type 2D report created.



Give the gift of DNA discovery.

Gift a kit

Refer friends, earn rewards.

Get reward

ANCESTRY

Ancestry Overview
All Ancestry Reports
Ancestry Composition
DNA Relatives
Order Your DNA Book

HEALTH & TRAITS

Health & Traits Overview
All Health & Traits Reports
My Health Action Plan
Health Predisposition
Pharmacogenetics
Carrier Status
Wellness
Traits

RESEARCH

Research Overview
Surveys and Studies
Edit Answers
Publications

FAMILY & FRIENDS

View all DNA Relatives
Family Tree
Your Connections
GrandTree
Advanced DNA Comparison