

Health > Health Predisposition

**HOME ANCESTRY**  **HEALTH & TRAITS** 

RESEARCH

**FAMILY & FRIENDS** 





Print

# **Alpha-1 Antitrypsin Deficiency**

AAT deficiency is a genetic condition that can lead to lung and liver disease. It is caused by decreased levels of the alpha-1 antitrypsin (AAT) protein. This test includes the two most common variants linked to this deficiency.

Overview

**Scientific Details** 

**Frequently Asked Questions** 

# deficiency. People with this result have an increased risk of developing lung and liver disease related to AAT deficiency.

Smoking, drinking excessive amounts of alcohol, and other factors can also affect your risk.

Jamie, you have two copies of a genetic variant linked to AAT

Variant detected



# in the SERPINA1 gene



### This test does not diagnose AAT deficiency or any other health conditions.

How To Use This Test

condition runs in your family, you think you might have this condition, or you have any concerns about your results.

Please talk to a healthcare professional if this

**Review the Genetic Health Risk tutorial See Scientific Details** 

**See Frequently Asked Questions** 

deficiency.

**Intended Uses** 

Limitations

Does not test for all possible variants linked to AAT deficiency.

Tests for the PI\*Z and PI\*S variants in the SERPINA1 gene linked to AAT

# **Ethnicity Considerations**

people of European descent.

• The variants included in this test are most common and best studied in

It is important to discuss this result with a healthcare professional.

You have an increased risk of developing lung and liver

disease related to AAT deficiency based on your genetic

result.

**See Scientific Details** 



We detected two copies of a variant called PI\*Z.

You inherited one copy of this variant from each of your parents.

influence your risk.

However, some people with this result do not develop lung or liver disease.


Smoking, drinking excessive amounts of alcohol, and other factors can also

will also develop liver disease.

children.

disease.






If your partner has a variant linked to AAT deficiency, your children could inherit

a combination of variants that may put them at risk of developing lung or liver

Consult with a healthcare professional before making any major lifestyle changes.

other factors.



Because it is a genetic condition, AAT deficiency is present

life, and age of onset is strongly affected by smoking. Some

people may never have symptoms of lung disease, especially

at birth. Symptoms of lung disease usually appear later in

if they don't smoke. Liver problems may develop anytime

Potential signs and symptoms

Shortness of breath and wheezing

• Lung disease, including emphysema

Liver disease, including cirrhosis

Recurrent lung infections

When it develops

from infancy to adulthood.

Chronic cough

**Smoking** 

People with AAT deficiency who smoke are

more likely to develop lung disease and to





**Smoking** 

Personal or family

history of lung







Occupational and

other exposures

**Certain infections** 

### AAT deficiency is most common in people of Northern European descent. In the U.S., 1 in 3,000-5,000 people has this condition.

How it's treated

How common is the condition?

symptoms, treatment focuses on management of lung and liver problems. Direct replacement of the AAT protein into the blood may be used to slow the progression of lung disease. Lung and liver transplants may be beneficial in

There is currently no known cure. People with AAT deficiency


are encouraged to avoid smoking, limit alcohol consumption,

and consider getting certain vaccinations. For those with

Read more at: National Heart, Lung, and Blood Institute GeneReviews MedlinePlus



It is important to discuss this result with a healthcare professional.




If you have a personal or family history of lung or liver disease, consult with a



**Learn more** 

healthcare professional.



**Membership Terms** 

**ANCESTRY** 

Earn up to \$20 for each friend you

All Ancestry Reports **Ancestry Composition DNA Relatives** 

**Ancestry Overview** 

Health & Traits Overview All Health & Traits Reports

**HEALTH & TRAITS** 

**Cookie Policy** 


Surveys and Studies **Edit Answers Publications** 

**RESEARCH** 

Research Overview

View all DNA Relatives Family Tree **Your Connections** GrandTree Advanced DNA Comparison

**FAMILY & FRIENDS** 



©2023 23andMe, Inc.

refer.

**Terms of Service** 


Order Your DNA Book


**Privacy Statement** 

My Health Action Plan Health Predisposition **Carrier Status** Wellness **Traits** 

Help

Resources





Health > Health Predisposition

Print

# Alpha-1 Antitrypsin Deficiency AAT deficiency is a genetic condition that can lead to lung and liver disease. It is caused by decreased levels

of the alpha-1 antitrypsin (AAT) protein. This test includes the two most common variants linked to this deficiency.

**Scientific Details Overview Frequently Asked Questions** 

AAT deficiency is caused by variants in the SERPINA1 gene.

## The SERPINA1 gene contains instructions for making a protein called alpha-1 **Chromosome 14**

**SERPINA1** 

# antitrypsin (AAT). This protein is made in the liver, but is transported to the

damage, and the liver can become damaged as well. Read more at MedlinePlus

Variants Detected

Genotype\*

lungs where it has a protective function. Certain variants in SERPINA1 result in

too little AAT protein getting transported to the lungs, and more AAT protein


getting trapped in the liver. As a result, the lungs are less protected from

Marker Tested

**Gene: SERPINA1** 

Marker: rs28929474

PI\*Z



View All Tested Markers

### **Biological explanation** Variant copy from one Variant copy from Typical vs. variant DNA sequence(s) of your parents your other parent

**Additional Information** 

	<ul> <li>Percent of 23andMe customers with variant</li> <li>References [ 2, 3, 5, 8, 11, 15, 16, 17, 19, 20, 23 ]   ClinVar</li> </ul>		
*This test cannot distinguish which capy you received from which parent. This test also	cannot determine whether multiple variants, if detected, were inherited from only one parent or from		
both parents. This may impact how these variants are passed down.	cannot determine whether multiple variants, if detected, were inherited from only one parent of from		
23andMe always reports genotypes based on the 'positive' strand of the human genome strand.	e reference sequence (build 37). Other sources sometimes report genotypes using the opposite		

descent who have the variants included in this test. Estimates for other ethnicities are not currently available. Keep in mind that other risk factors — including smoking, drinking excessive amounts of alcohol, and having nonalcoholic fatty liver disease (NAFLD) — can increase the risk of developing lung and severe liver disease, regardless of genetics.

Average serum

AAT levels, µM/L

(5th to 95th %ile)

Test Interpretation

This report provides information about the risk of developing lung and liver disease in people of European

### genotype and a health condition. PI\*MS 33 (18-52) For certain genotypes, quantitative risk

estimates may not be available. including emphysema, Variants in the SERPINA1 gene can affect AAT due to AAT deficiency. protein levels differently. Severe AAT deficiency is defined by AAT levels below 11 µM/L. Lung diseases such as emphysema and COPD are most commonly associated with AAT levels below this protective threshold. This

Genotype (1)

This is not a complete list of other factors.

The factors described here include the most

associated with lung or liver disease in people

with AAT deficiency. Other factors not listed

here may also influence risk for lung or liver

Consult with a healthcare professional before

disease in people with the condition.

making any major lifestyle changes.

common and well-established risk factors

**Health Risk Estimates** 

that identify an association between a

Risk estimates are based on clinical studies

table provides AAT protein levels associated

only, and does not indicate a person's actual

protein levels.

**References** [ 1, 16, 20 ]

with each genotype for informational purposes

Consider talking to a healthcare professional if

you have any concerns about your results.

/ AAT levels below 11					
ich as emphysema and only associated with AAT tive threshold. This tein levels associated informational purposes tate a person's actual	PI*SS	28 (20–48)	Not likely at increased risk of developing COPD, including emphysema, due to AAT deficiency.	Not likely at increased risk of developing severe liver disease, including cirrhosis, due to AAT deficiency.	
ealthcare professional if about your results.	PI*MZ	25.4 (15–42)	Not likely at increased risk of developing COPD, including emphysema, due to AAT deficiency. However, smokers with this genotype have an increased risk.	Not likely at increased risk of developing severe liver disease, including cirrhosis, due to AAT deficiency. However, excessive alcohol consumption and having nonalcoholic fatty liver disease (NAFLD) can increase risk.	
	PI*SZ	16.5 (10-23)	Not likely at increased risk of developing COPD, including emphysema, due to AAT deficiency.  However, scientists estimate that 20–50% of smokers with this genotype will develop signs of emphysema during their lifetime.	Not likely at increased risk of developing severe liver disease, including cirrhosis, due to AAT deficiency. However, excessive alcohol consumption and having nonalcoholic fatty liver disease (NAFLD) may increase risk.	
	PI*ZZ	5.3 (3.4–7.0)	Increased risk of developing COPD, including emphysema, due to AAT deficiency. Scientists estimate that greater than 80% of people with this genotype will develop signs of emphysema during their lifetime.	Increased risk of developing severe liver disease due to AAT deficiency. Scientists estimate that people with this genotype have a 30–40% chance of developing cirrhosis after the age of 50.	
			gher risk of developing lung and er factors.		
ict of other feeters	Other Factors			References	
ist of other factors. ere include the most	Smoking			[ 1, 12, 20 ]	
liver disease in people her factors not listed e risk for lung or liver	People with at least one copy of the PI*Z <u>variant</u> are more likely to develop lung disease if they smoke. People with AAT deficiency who smoke typically start to experience the symptoms of lung disease between 40 and				

Risk estimates for developing lung and liver disease in people of European descent

**Lung Disease** 

Not likely at increased

risk of developing COPD,

**Liver Disease** 

Not likely at increased

risk of developing severe

liver disease, including

cirrhosis, due to AAT

deficiency.

## to industrial gases, metal fumes, and mineral dust may lead to a faster decline in lung function in people with AAT deficiency. Exposure to pollutants from kerosene heaters on a regular basis may also increase the chances of developing lung disease related to AAT deficiency. The effects

Personal or family history of lung disease

develop the symptoms of lung disease.

Occupational and other exposures

People with AAT deficiency whose siblings suffer from lung disease are more likely to develop lung disease themselves. This may be due to genetic and/or environmental factors. **Certain infections** Diseases like the flu can damage the lungs, and diseases like hepatitis A and B can damage the liver. Yearly immunization against influenza (a virus that causes the flu) and immunization against pneumococcus (a bacterium

that causes respiratory infections) are generally recommended for people

with AAT deficiency. This can prevent lung disease from getting worse.

Immunizations against the viruses hepatitis A and B, which cause liver

People with AAT deficiency who have a personal history of lung problems

later in life. The risk of lung disease can also depend on family history.

such as asthma or wheezing are more likely to develop severe lung disease

50 years of age. In contrast, non-smokers with AAT deficiency may not

experience symptoms until their 60s, and some non-smokers will never

A small number of research studies — mostly looking at men working in

of occupational and other exposures on symptoms of lung disease in

people with AAT deficiency are still not fully understood.

construction or farming — suggests that prolonged occupational exposure

(NAFLD). This means reducing risk for NAFLD may help lower the chances of developing cirrhosis. Factors like maintaining a healthy weight and keeping blood sugar and cholesterol levels in the healthy range can help reduce the risk for NAFLD.

**Certain health conditions** 

**Excessive alcohol consumption** 

The 23andMe PGS Genetic Health Risk Report for Alpha-1 Antitrypsin Deficiency is indicated for reporting of the PI\*Z and PI\*S variants in the SERPINA1 gene. This report describes if a person has variants associated with AAT deficiency and a higher risk for lung or liver disease, but it does not describe a person's overall risk of developing lung or liver disease. This report is most relevant for

Test Details

## Accuracy was determined by comparing results from this test with results from sequencing. Greater than 99% of test results were correct. While unlikely, this test may provide false positive or false negative results. For more details on the analytical performance of this test, refer to the package insert.

Pulmon Dis. 12:1683-1694.

Obstr Pulm Dis. 7(3):260-271.

**Test Performance Summary** 

**Indications for Use** 

people of European descent.

**Special Considerations** 

**Clinical Performance** 

**Analytical Performance** 

SERPINA1 gene.

Dis. 12:561-569. \ 4. Chalasani N et al. (2018). "The diagnosis and management of nonalcoholic fatty liver disease: Practice guidance from the American Association for the Study of Liver Diseases." Hepatology. 67(1):328-357.

Testing for genetic variants associated with AAT deficiency is recommended under certain

More than 95% of all cases of AAT deficiency are caused by the PI\*Z and PI\*S variants in the

circumstances by several health professional organizations, including the American Thoracic Society.

References

1. American Thoracic Society. et al. (2003). "American Thoracic Society/European Respiratory Society statement: standards for the diagnosis and

2. Blanco I et al. (2017). "Alpha-1 antitrypsin Pi\*SZ genotype: estimated prevalence and number of SZ subjects worldwide." Int J Chron Obstruct

3. Blanco I et al. (2017). "Alpha-1 antitrypsin Pi\*Z gene frequency and Pi\*ZZ genotype numbers worldwide: an update." Int J Chron Obstruct Pulmon

8. Hersh CP et al. (2004). "Chronic obstructive pulmonary disease in alpha1-antitrypsin PI MZ heterozygotes: a meta-analysis." Thorax. 59(10):843-9.

management of individuals with alpha-1 antitrypsin deficiency." Am J Respir Crit Care Med. 168(7):818-900.

5. Dahl M et al. (2005). "The protease inhibitor PI\*S allele and COPD: a meta-analysis." Eur Respir J. 26(1):67-76.

[1]

- 6. Fregonese L et al. (2008). "Hereditary alpha-1-antitrypsin deficiency and its clinical consequences." Orphanet J Rare Dis. 3:16. 7. Hamesch K et al. (2020). "Non-Invasive Assessment and Management of Liver Involvement in Adults With Alpha-1 Antitrypsin Deficiency." Chronic
- 9. Köhnlein T et al. (2010). "Diagnostic delay and clinical modifiers in alpha-1 antitrypsin deficiency." Ther Adv Respir Dis. 4(5):279-87.
- Change Log

10. Lieberman J et al. (1986). "Alpha 1-antitrypsin Pi-types in 965 COPD patients." Chest. 89(3):370-3.

Information about liver disease risk was updated for people with certain Nov. 3, 2021 genotypes.

April 17, 2017	Alpha-1 Antitrypsin Deficiency report created.				
	ANCESTRY	HEALTH & TRAITS	RESEARCH	FAMILY & FRIENDS	
	Ancestry Overview	Health & Traits Overview	Research Overview	View all DNA Relatives	

refer. **Get \$20** 

23andMe

Earn up to \$20 for each friend you

**Date** 

**DNA Relatives** Order Your DNA Book

Change

**Carrier Status** Wellness

GrandTree

Advanced DNA Comparison

View all DNA Relatives Family Tree **Your Connections** 

disease, may also be recommended by a healthcare professional. [ **7**, **21** ] People with at least one copy of the PI\*Z variant are more likely to develop severe liver disease if they drink excessive amounts of alcohol. [4, 7, 21] People with at least one copy of the PI\*Z variant are more likely to develop severe liver disease if they also have nonalcoholic fatty liver disease

**Warnings and Limitations** 

• This test does not cover all variants that

• If you are concerned about your results,

consult with a healthcare professional.

See the **Package Insert** for more details on

\* Variants not included in this test may be very rare,

may not be available on our genotyping platform, or

use and performance of this test.

may not pass our testing standards.

could cause this condition.\*

conditions.

[ **1**, **18** ]

[1]

[ **1**, **9** ]

• This test does not diagnose any health • Share results with your healthcare professional for any medical purposes.

See all references >

**Traits** 



**HOME** 









Health > Health Predisposition



## Alpha-1 Antitrypsin Deficiency

AAT deficiency is a genetic condition that can lead to lung and liver disease. It is caused by decreased levels of the alpha-1 antitrypsin (AAT) protein. This test includes the two most common variants linked to this deficiency.

> **Scientific Details** Overview

**Frequently Asked Questions** 

## Alpha-1 Antitrypsin Deficiency

What does this test do?	~
What does this test <b>not</b> do?	~
The report says the variants included in this test are most common and best studied in people of <b>European</b> descent. What if I'm not of European descent?	~
Where can I learn more about alpha-1 antitrypsin deficiency, support groups, and other resources?	~
My report says <b>two copies of a variant</b> called <b>PI*Z</b> were detected. What does this mean?	~
What does increased risk mean?	~
My report says <b>two copies of a variant</b> called <b>PI*Z</b> were detected. What are some things I could do?	~
How could my result affect my family?	~

Have more questions? Check out our Customer Care Help Center.



Earn up to \$20 for each friend you refer.

Get \$20

**ANCESTRY** 

**DNA Relatives** 

**Ancestry Overview** All Ancestry Reports

**Ancestry Composition** 

Order Your DNA Book

**HEALTH & TRAITS** 

Health & Traits Overview All Health & Traits Reports

Health Predisposition

My Health Action Plan

**Carrier Status** Wellness

Traits

**RESEARCH** 

Research Overview Surveys and Studies

**Edit Answers** 

**Publications** 

**FAMILY & FRIENDS** 

View all DNA Relatives

Family Tree

**Your Connections** 

GrandTree

Advanced DNA Comparison





